## **ANNEXES**

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**(NOTE D'APPLICATION PERKIN ELMER)**

# Uptake Rates for Perkin-Elmer Diffusive Tubes

## Application Note

This information has been reproduced from Issue 7 of The Diffusive Monitor with the kind permission of the editor, Dr. R. H. Brown. "The Diffusive Monitor" is a free-issue magazine/newsletter focusing on "current news in the field of diffusive sampling in the workplace." It is issued by Working Group 5 of the UK Health and Safety Executive (HSE/CAR) and presents many articles and points of interest to the ATD user.

### Diffusive Uptake Rates on Perkin-Elmer Sorbent Tubes

NOTE: This list has been compiled from sources available to CAR/Working Group 5. It is not

complete, and although every care has been taken in its preparation, no guarantee can be given of its accuracy. Nor does the list imply a preference for this type of sampler.

Sampling rates are given for the standard Perkin-Elmer tube without the membrane, except where noted.





Separate lists are given for workplace (approx 8 hours) and environmental (approx 4 weeks) exposure periods and exposure concentrations typical of these applications.

Users of the data are strongly advised to consult the original source material to determine the level of confidence and range of applicability of the values.

### Workplace Applications

#### Hydrocarbons

Compound	Sorbent	Level	Uptake Rate**	Source
1,3-butadiene	Mol. Sieve 13X	A	1.3	14
n-pentane	Chrom. 106	A	1.46	15
	Carbopack B*	B	1.77	2
n-hexane benzene	Chrom. 106	A	1.77	15
	Tenax TA*	A	1.3	16
	Porapak Q	A	1.37	16
	Tenax GR	B	1.81	17
n-heptane	Chrom. 106*	B	1.72	4
	Chrom. 106	A	1.95	15
	Tenax TA*	A	1.77	15
toluene	Carbopack B	B	1.94	2
	Tenax TA*	B	1.67	4
	Tenax GR	B	2.12	17
	Chrom. 106	B	1.94	4
n-octane	Chrom. 106	B	2.06	2
	Chrom. 106	A	2.13	15
	Tenax TA*	A	2.00	15
xylene	Tenax TA*	B	1.82	4
	Chrom. 106	B	2.10	4
	Tenax GR	B	2.48	17

**PERKIN ELMER**

*Hydrocarbons, continued*

Compound	Sorbent	Level	Uptake Rate**	Source
ethyl benzene	Tenax TA*	B	2.0	4
	Tenax GR	B	2.43	17
	Chrom. 106	B	1.9	4
	Porapak Q	D	2.38	18
styrene	Tenax TA*	A	2.0	19
	Chrom. 106	B	2.15	4
n-nonane	Chrom. 106	A	2.40	15
	Tenax TA*	A	2.12	15
iso-propyl benzene	Chrom. 106	C	2.26	15
	Tenax TA*	C	2.26	15
trimethyl-benzene	Chrom. 106	C	2.37	15
	Tenax TA*	C	2.37	15
n-decane	Tenax TA	A	2.3	5
cumene	Porapak Q	D	2.5	18

*Halogenated Hydrocarbons*

Compound	Sorbent	Level	Uptake Rate**	Source
methyl chloride	Spherocarb	B	1.3	18
vinyl chloride	Spherocarb	B	2.0	18
1,1-dichloro-ethene	Spherocarb	B	2.5	18
trichloro trifluoro ethane	Chrom. 102	B	3.5	18
chloro trifluoro methane	Chrom. 102	B	1.8	18
dichloro methane	Chrom. 106	B	1.5	4
	Chrom. 102	B	1.56	18
1,2-dichloro ethane	Chrom. 102	B	1.9	18
fluothane	Chrom. 102	B	3.6	18
halothane	Tenax TA	B	2.59	6
enflurane	Tenax TA	B	2.29	6
isoflurane	Tenax TA	B	2.20	6
bromoethane	Chrom. 106	E	2.45	20
trichloro-methane (chloroform)	Tenax GR	B	2.18	17
	Chrom. 102	B	2.35	18
tetrachloro-methane (carbon tetrachloride)	Tenax GR	B	3.72	17
	Chrom. 102	B	2.87	18
trichloro ethene	Chrom. 106	B	2.5	4
	Chrom. 102	B	2.3	18
1,1,1-trichloro-ethane	Chrom. 106	B	2.3	4
	Chrom. 102	B	2.3	18
	Tenax GR	B	2.92	17
tetrachloro ethene	Chrom. 106	B	3.1	4
	Tenax TA	B	2.8	4
	Chrom. 102	B	2.6	18
epichlorohydrin	Chrom. 106	E	2.45	20

*Halogenated Hydrocarbons, continued*

Compound	Sorbent	Level	Uptake Rate**	Source
perfluoro-dimethyl cyclobutane***	Carbotrap	B	15 ml./hour	7
perfluoro-methyl-cyclopentane***	Carbotrap	B	15 ml./hour	7
perfluoro-methyl-cyclohexane***	Carbotrap	B	15 ml./hour	7

*Esters and Glycol Ethers*

Compound	Sorbent	Level	Uptake Rate**	Source
ethyl acetate	Chrom. 106*	B	1.77	4
	Tenax TA	B	1.43	4
n-butyl acetate	Tenax TA	B	2.26 (8 hour)	8
methyl methacrylate	Porapak Q	B	2.0	18
butyl acrylate	Tenax TA	B	2.6 (8 hour)	8
2-methoxy-ethanol	Porapak Q	A	1.5	21
2-ethoxy-ethanol	Tenax	A	1.8	21
2-methoxyethyl acetate	Porapak Q	A	2.8	21
2-ethoxyethyl acetate	Chrom. 106	B	2.3	4
	Tenax TA	B	2.1	4
butoxyethanol	Chrom. 106	B	2.1	4
	Tenax TA	B	1.9	4
methoxy-propanol	Chrom. 106*	B	1.85	4
	Tenax TA	B	1.52	4
Methyl oxitol	Chrom. 106	B	2.1	9
2-butoxyethyl acetate	Tenax	A	2.8	28

*Aldehydes and Ketones*

Compound	Sorbent	Level	Uptake Rate**	Source
methyl isobutyl ketone (2-methyl-4-pentanone)	Tenax TA	B	1.71	4
	Chrom. 106*	B	2.01	4
Cyclohexanone	Tenax TA	D	2.3	5
furfural	Tenax TA	A	2.5	10

*Alcohols*

Compound	Sorbent	Level	Uptake Rate**	Source
propan-2-ol (isopropanol)	Spherocarb	C	2.0	22

Miscellaneous

Compound	Sorbent	Level	Uptake Rate**	Source
acrylonitrile	Porapak N	A	1.35	23
acetonitrile	Porapak N	A	1.0 (2 hrs.)	11
	Porapak N	A	0.8 (8 hrs.)	11
propionitrile	Porapak N	A	1.4 (2 hrs.)	11
	Porapak N	A	1.3 (8 hrs.)	11
carbon disulphide	Spherocarb	A	2.6	5
nitrous oxide	Mol. Sieve 5A <sup>¶¶</sup>	B	1.25	12
ethylene oxide	Spherocarb	B	1.6	13
1,4-dioxane	Spherocarb	C	3.0	22

Environmental Applications

Hydrocarbons

Compound	Sorbent	Level	Uptake Rate**	Source
benzene	Tenax TA	B	0.88	24
	Chrom. 106	B	1.89	24
	Tenax GR	B	0.8	1
toluene	Tenax TA	B	1.2	3
xylene	Tenax TA	B	1.91	3
trimethylbenzene	Tenax TA	B	2.65	3
decane	Tenax TA	B	2.96	3
undecane	Tenax TA	B	3.38	3

Notes to the Tables

Footnotes

\* means preferred sorbent

\*\* units are ng.ppm<sup>-1</sup>.min<sup>-1</sup>; most values are for 4-8 hrs: environmental values are for 4 weeks.

\*\*\* A nickel disk, rather than the conventional stainless steel gauze, was used to support the Carbotrap sorbent material during method validation for these perfluorocarbon tracer gases. The uptake rates may not be applicable to samples using conventional steel gauzes.

¶ sorbent not recommended because of varying uptake rate

¶¶ rate varies predictably with exposure dose

Levels

level A = validation equivalent to CEN level 1A  
 level B = partial validation - see EN 482  
 level C = calculated - ideal value  
 level D = calculated from dynamic breakthrough volume  
 level E = calculated from sorption isotherm

Sources/References

1. Peters, R and Hafkenschied, Th. *The Diffusive Monitor*, 7, 8-9 (1995)

2. Van den Hoed N., van Asselen O.L.J. A computer model for calculating effective uptake rates of tube-type diffusive air samplers. *Ann. Occup. Hyg.*, 35, 273-285. (1991)

3. Brown V.M., Crump D.R., Gardiner D. Yu C.W.F. Long term diffusive sampling of volatile organic compounds in indoor air. *Environmental Technology*, 14, 771-777. (1993).

4. Wright M.D. Diffusive uptake rates for the Perkin-Elmer tube - BCR air sampling intercomparisons at Vito (Mol. Belgium) February 1991 - April 1992. UK Health and Safety Executive Report Ref. IR/U 1A/93/3.

5. UK Health and Safety Executive, Occupational Medicine and Hygiene Laboratory. Diffusion methods for the determination of gas and vapour pollutants in air. Internal Health and Safety Executive Guidance Document, 1988.

6. Gray W.M., O'Sullivan J., Houldsworth H.B., Musgrove, N.T. The use of diffusive samplers to monitor occupational exposure to waste anaesthetic gases. CEC Publ No. 10555EN, Brussels, Luxembourg (1987).

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8. Rosmanith P. Determination of some uptake rates for Perkin-Elmer ATD 50 diffusive sampling: Diffusive sampling-an alternative approach to workplace air monitoring. CEC Publ. No. 10555E.

9. Glover J.H. Thermal desorption in industrial hygiene and environmental analysis. Spantech House, South Godstone, Surrey, England, Spantech Publishers, 1991.

10. Patel S. Method validation for atmospheric monitoring - an illustration with furfural. 1988 Thermal Desorption Symposium, Windsor UK.

11. Wood M. The use of the Perkin-Elmer passive sampler and ATD 50 automatic thermal desorber in the measurement of atmospheric concentrations of organic nitriles. *Ann. Occup. Hyg.*, 29, 399-413 (1985)

- |  |   |   |
|--|---|---|
| <p>12. Cox P.C., Brown R.H. A personal sampling method for the determination of nitrous oxide exposure. <i>Am. Ind. Hyg. Assoc. J.</i>, 45, 345 (1984)</p> <p>13. HSE Unpublished data.</p> <p>14. MDHS 63</p> <p>15. MDHS 66</p> <p>16. MDHS 50</p> <p>17. TNO, Netherlands (contact: Dr T. Hafkenscheid)</p> | <p>18. ICI, UK (contact: Mr. S Bennett)</p> <p>19. MDHS 43</p> <p>20. Unwin, J. Determination of diffusive uptake rates for tube-type samplers using sorption isotherms and computer modelling. <i>The Diffusive Monitor</i>, 6.6 (1993)</p> <p>21. BP International, UK (contact: Dr K J Saunders)</p> <p>22. Rank Xerox, UK</p> | <p>23. MDHS 55</p> <p>24. Saunders, KJ The use of diffusive samplers for measuring benzene at ppb levels in the environment. <i>The Diffusive Monitor</i>, 6.4-5 (1993)</p> <p><i>Note: HSE Internal Reports may be obtained from the HSE Information Centre, Broad Lane, Sheffield, UK. Tel +44 1142 892345, FAX +44 1142 892333</i></p> |
|--|---|---|

## PERKIN ELMER

### The Perkin-Elmer Corporation

761 Main Avenue, Norwalk, CT 06859-0012 U S A  
 Tel: (800) 762-4000 or (203) 762-4000 Fax: (203) 762-4222  
 Internet: <http://www.perkin-elmer.com>

All analytical instruments and systems manufactured by Perkin-Elmer are developed and produced under the quality requirements of ISO 9001. Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation.

Printed in U S A © 1996 The Perkin-Elmer Corporation  
 Order No GCA-72 March 1996 KG039602  
 Replaces TDD-2